

HAEQUEO	Product Technical Dossier	Rev.02	
	cross-linked sodium hyaluronate 2,5% intra-dermal use	Issued on: 30.01.2016	

Sterilisation process

At the end of the productive process the filling of syringes is carried out and therefore also the steam sterilisation in autoclave.

The sterilisation process is performed in accordance with the harmonized standards and rules EN 554 and ISO 17665.

The relative registrations to parameters for the sterilisation process used during each process are stored for 5 years after the LOT expiry date to which they refer to. The relative data for the sterilisation process validation is stored for 10 years.